

A. D. Darbal<sup>1,4</sup>, R. D. Narayan<sup>1</sup>, C. Vartuli<sup>2</sup>, G. Lian<sup>2</sup>, R. Graham<sup>3</sup>, F. Shaapur<sup>3</sup>, S. Nicolopoulos<sup>4</sup> and J. K. Weiss<sup>1</sup>

<sup>1</sup> AppFive LLC, Tempe, Arizona, USA. <sup>2</sup> Texas Instruments, Dallas, Texas, USA <sup>3</sup> Nanotem Inc, Scottsdale, Arizona, USA <sup>4</sup> NanoMEGAS SPRL, Brussels, Belgium

## INTRODUCTION

- Measurement of strain with high spatial resolution and high precision in semiconductor devices is critical to monitor the designed and unintended strain distributions.
- Use of spot diffraction patterns with nanobeam illumination gives higher spatial resolution than other TEM techniques[1].
  - Experiment is relatively simple.
- Technique is made possible by beam precession
  - Improves quality of diffraction patterns..

[1] D Cooper *et al.*, Journal of Physics: Conference Series **326** 012025 (2011).

## LIMITATIONS OF EXISTING TEM STRAIN MEASUREMENT METHODS

Technique	Advantages	Limitations
Convergent beam electron diffraction	High spatial, strain sensitivity	<ul style="list-style-type: none"> <li>• Needs the sample to be relatively thick (&gt;150 nm)</li> <li>• Sample needs to be oriented away from a low index axis.</li> <li>• Very sensitive to strain relaxation</li> </ul>
Dark field holography	High spatial resolution (5 nm), large field of view (1 $\mu\text{m} \times 1 \mu\text{m}$ )	<ul style="list-style-type: none"> <li>• Requires unstrained reference with identical crystallographic orientation area close to strained region</li> </ul>
High resolution imaging	High spatial resolution (< 1 nm)	<ul style="list-style-type: none"> <li>• Limited field of view (100 <math>\text{nm}^2 \times 100 \text{nm}^2</math>)</li> <li>• Stringent requirements on specimen quality</li> </ul>


## CONVENTIONAL NANOBEAM DIFFRACTION

- Acquire spot diffraction patterns from strained and unstrained regions using a quasi-parallel nanoprobe (<5 nm)
- Use measured shift in spot positions to calculate strain
- Experiment is relatively straightforward


### Limitations

- Presence of strong dynamical effects lead to rapid changes in spot intensities with small thickness and orientation changes
- Strong dependence of spot intensities on changes in local thickness and orientation makes automated analysis challenging
  - Requires manual intervention in identifying spot positions
- Inadequate sampling of higher order reflections limits the accuracy

## PRECESSION ELECTRON DIFFRACTION



Portillo, J., *et al.* (2010). *Materials Science Forum* **23**, 1-7




Diffraction patterns from two points 120 nm apart from Si/SiGe multilayered specimen

- No particular beam is strongly diffracted – reducing strong dynamical effects
- Insensitive to small thickness and orientation changes
- Number of spots increases – better sampling of higher order spots

## STRAIN MEASUREMENT ANALYSIS


- Diffraction patterns from strained region are matched against a reference pattern.
- Reference pattern from unstrained region.
- Correlation distance used as the metric for fitting reference to strained patterns.
- Results include strain in x and y-directions and shear (not shown).
  - Relative to x-direction specified by user.



## RESULTS


- Data acquired with Zeiss Libra L200 TEM.
  - Field Emission Gun (FEG)
  - Scanning TEM (STEM) mode.
  - NanoMEGAS DigiSTAR unit for precession and descanning of beam.
- Positive percentage strain values correspond to tensile strain, negative values compressive.
- Precision of strain measurement is 0.02% in profile below.
  - Precision can depend on characteristics of data, such as resolution of spot patterns.

### 1. Strain profile of an Si/SiGe layer.




- Strain in x-direction is near zero, indicating a coherent interface.

### 2. Strain maps from a AlGaIn/GaN HEMT



- For such devices, tensile strain is expected to be asymmetric on different sides of the gate, and this is seen here.
- Compressive strain perpendicular to the interface (x-direction) of ~1% is seen in the AlGaIn region.

### 3. Strain maps from the Si region of a pMOS device.



- x and y-directions aligned with [220] and [002] directions in Si.
- Localized biaxial tensile strain close to contact edges.